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Vertrieb/ Distribution:

SERAG-WIESSNER GmbH & Co. KG
95119 Naila / Deutschland

Sterile Polydioxanone Suture With Needle
Intended use

This device is intended to fixate sub dermal tissue in an elevated position in plastic and reconstructive surgery. Polydioxanone suture provides wound support for a longer period compared to other synthetic absorbable sutures. It also prolongs the stimulation time that increases the treatment effect on the body.

Device characteristics

´Sterile Polydioxanone Suture with Needle´ consists of needle, needle cap, hub, sponge and polydioxanone (PDO) suture and is sterilized by EO gas. The safety and effectiveness of polydioxanone (PDO) is proven.

Things to do before use

- Before using the product, the condition of the patient must be stable
- Prepare disinfection of the area to be treated, e.g. with alcohol cottons.
- Prepare the single use ´Sterile polydioxanone suture with needle´ that is packed sterile.
- Check whether the model is correct.
- Check the minimum durability date.
- Check whether the package is damaged.
- Use only after understanding completely how to use this product.

How to use the product

- Check the content and condition of the packaging.
- Get to know the instructions for use thoroughly.
- Use the products only if they´re not damaged and packed sterile.
- Remove the package of the product, and carefully take out ´Sterile polydioxanone suture with needle´.
- Confirm the integrity of the product.
- Thoroughly disinfect the areas to be treated with ´Sterile polydioxanone suture with needle´.
- Insert the ´Sterile polydioxanone suture with needle´ under the skin or muscle to be treated.
- Slowly and gently remove the needle or cannula and leave the thread in the appropriate skin layer.

How to manage after use

- Discard the product if it was opened.
- This product is a single use, must not be reused.
- Discard the used product into an appointed safe place.

Warning

Patients on whom the sutures should not be used:

- Patients with malnutrition.
- Debilitated patient.
- Patients with hypertrophic scar formation.
- Patients that are treated with anticoagulants or thrombocyte aggregation inhibitors (e.g.ASS).
- Patients with infected or inflamed skin areas.
- Patients with autoimmune diseases who are getting an immunotherapy.
- Patients with known sensitivity to polydioxanone.
- Pregnant and breastfeeding women.
- Patients under 18 years.

Side effects

Generally, the PDO sutures get tolerated very well. Very rarely there can occur injection induced skin reactions such as light pain, tensioned and pressured feeling, stinging, itching, hematoma or light swellings. This reactions ease after a short period of time.

Precautions

- This product should be used by suitably trained practitioners.
- Do not use the product if the minimum durability date has exceeded.
- Do not use the product if the packaging is damaged.
- This product must not be used for other purposes than the intended one.

If the area to be treated is inflamed or infected, the product should not be used.

Preservation method

Store the product in a clean and dry room. Avoid direct sunlight and high moisture.

Contraindication

If the patient suffers from cardiac disorder, high blood pressure, diabetes, panic disorder, depression, etc., there has to be further discussion with the doctor.

Prohibition

- Before the operation, for about a week, stop taking anticoagulants (e.g. Aspirin, Marcumar etc.)
- Before and after the operation, stop smoking and drinking.
- Unforeseen diseases such as a cold or fever may cause problems regarding anesthesia; the treating doctor has to be informed about such diseases.

Expiration date

2 years from the date of sterilization.

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